

For Research Use Only

# TBC1D17 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60975-2-PBS



## Basic Information

<b>Catalog Number:</b> 60975-2-PBS	<b>GenBank Accession Number:</b> BC003516	<b>Purification Method:</b> Protein A Magarose purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 79735	<b>CloneNo.:</b> 3F5C2
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9HA65	
<b>Isotype:</b> IgG3	<b>Full Name:</b> TBC1 domain family, member 17	
<b>Immunogen Catalog Number:</b> AG14358	<b>Calculated MW:</b> 648 aa, 73 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA, Sample test

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sample test of MP51406-1, TBC1D17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60975-1-PBS. Detection antibody: 60975-2-PBS.



Cytometric bead array standard curve of MP51406-1, TBC1D17 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60975-1-PBS. Detection antibody: 60975-2-PBS. Standard: Ag14358. Range: 1.563-200 ng/mL