

For Research Use Only

LYVE1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60976-1-PBS



Basic Information

Catalog Number:

60976-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG0570

GenBank Accession Number:

BC026231

GeneID (NCBI):

10894

UNIPROT ID:

Q9Y5Y7

Full Name:

lymphatic vessel endothelial
hyaluronan receptor 1

Calculated MW:

322 aa, 35 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

4E4D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP51407-1, LYVE1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60976-1-PBS. Detection antibody: 60976-2-PBS.



Cytometric bead array standard curve of MP51407-1, LYVE1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60976-1-PBS. Detection antibody: 60976-2-PBS. Standard:Ag0570. Range: 0.781-200 ng/mL