

For Research Use Only

SLC7A1/CAT-1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60978-1-PBS



Basic Information

Catalog Number:

60978-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG5398

GenBank Accession Number:

BC063303

GeneID (NCBI):

6541

UNIPROT ID:

P30825

Full Name:

solute carrier family 7 (cationic amino acid transporter, y+ system), member 1

Calculated MW:

68 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2H4F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

