

For Research Use Only

SLC7A1/CAT-1 Monoclonal antibody, PBS Only (Detector)

Catalog Number:60978-2-PBS



Basic Information

Catalog Number: 60978-2-PBS	GenBank Accession Number: BC063303	Purification Method: Protein A Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 6541	CloneNo.: 2E11G4
Source: Mouse	UNIPROT ID: P30825	
Isotype: IgG2a	Full Name: solute carrier family 7 (cationic amino acid transporter, y+ system), member 1	
Immunogen Catalog Number: AG5398	Calculated MW: 68 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

