

For Research Use Only

SLC7A6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60980-2-PBS



Basic Information

Catalog Number:

60980-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4924

GenBank Accession Number:

BC028216

GeneID (NCBI):

9057

UNIPROT ID:

Q92536

Full Name:

solute carrier family 7 (cationic amino acid transporter, y+ system), member 6

Calculated MW:

515 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D11F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

