

For Research Use Only

# VSTM2L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60981-2-PBS



## Basic Information

Catalog Number:

60981-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22362

GenBank Accession Number:

BC033818

GeneID (NCBI):

128434

UNIPROT ID:

Q96N03

Full Name:

V-set and transmembrane domain  
containing 2 like

Calculated MW:

204 aa, 22 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B6H3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

