

For Research Use Only

# TAF13 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60984-1-PBS



## Basic Information

<b>Catalog Number:</b> 60984-1-PBS	<b>GenBank Accession Number:</b> BC121180	<b>Purification Method:</b> Protein A Magarose purification
<b>Concentration:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 6884	<b>CloneNo.:</b> 2G4C12
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q15543	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa	
<b>Immunogen Catalog Number:</b> AG18502	<b>Calculated MW:</b> 124 aa, 14 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:  
T: 4006900926      E: Proteintech-CN@ptglab.com      W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

