

For Research Use Only

TAF13 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60984-2-PBS



Basic Information

Catalog Number: 60984-2-PBS	GenBank Accession Number: BC121180	Purification Method: Protein A Magarose purification
Concentration: 1 mg/ml	GeneID (NCBI): 6884	CloneNo.: 2A1D8
Source: Mouse	UNIPROT ID: Q15543	
Isotype: IgG2a	Full Name: TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa	
Immunogen Catalog Number: AG18502	Calculated MW: 124 aa, 14 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

