

For Research Use Only

NMNAT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number:60994-1-PBS



Basic Information

Catalog Number: 60994-1-PBS	GenBank Accession Number: BC014943	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 64802	CloneNo.: 1B4G9
Isotype: IgG1	UNIPROT ID: Q9HAN9	
Immunogen Catalog Number: AG29593	Full Name: nicotinamide nucleotide adenylyltransferase 1	
	Calculated MW: 33 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51431-1, NMNAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60994-1-PBS. Detection antibody: 60994-2-PBS.



Cytometric bead array standard curve of MP51431-1, NMNAT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60994-1-PBS. Detection antibody: 60994-2-PBS. Standard: Ag29593. Range: 0.391-100 ng/mL.