

For Research Use Only

Connexin 32 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60996-2-PBS



Basic Information

Catalog Number: 60996-2-PBS	GenBank Accession Number: BC002805	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 2705	CloneNo.: 1F2C2
Isotype: IgG2a	UNIPROT ID: P08034	
Immunogen Catalog Number: AG0659	Full Name: gap junction protein, beta 1, 32kDa	
	Calculated MW: 32 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51433-1, Connexin 32 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60996-1-PBS. Detection antibody: 60996-2-PBS.



Cytometric bead array standard curve of MP51433-1, Connexin 32 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60996-1-PBS. Detection antibody: 60996-2-PBS. Standard: Ag0659. Range: 0.098-100 ng/mL.