

For Research Use Only

APOL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60998-2-PBS



Basic Information

Catalog Number:
60998-2-PBS

Source:
Mouse

Isotype:
IgG2b

Immunogen Catalog Number:
AG7622

GenBank Accession Number:
BC004395

GeneID (NCBI):
23780

UNIPROT ID:
Q9BQE5

Full Name:
apolipoprotein L, 2

Calculated MW:
37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

4G4D2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

