

For Research Use Only

APOL2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60998-3-PBS



Basic Information

Catalog Number:

60998-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7622

GenBank Accession Number:

BC004395

GeneID (NCBI):

23780

UNIPROT ID:

Q9BQE5

Full Name:

apolipoprotein L, 2

Calculated MW:

37 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G5G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51435-2, APOL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60998-3-PBS. Detection antibody: 60998-4-PBS.

Cytometric bead array standard curve of MP51435-2, APOL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60998-3-PBS. Detection antibody: 60998-4-PBS. Standard: Ag7622. Range: 0.781-100 ng/mL