

For Research Use Only

APOL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 60998-4-PBS



Basic Information

Catalog Number: 60998-4-PBS	GenBank Accession Number: BC004395	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 23780	CloneNo.: 4D6H2
Isotype: IgG1	UNIPROT ID: Q9BQE5	
Immunogen Catalog Number: AG7622	Full Name: apolipoprotein L, 2	
	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51435-2, APOL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60998-3-PBS. Detection antibody: 60998-4-PBS.

Cytometric bead array standard curve of MP51435-2, APOL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60998-3-PBS. Detection antibody: 60998-4-PBS. Standard: Ag7622. Range: 0.781-100 ng/mL