

For Research Use Only

DHX15 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 60999-1-PBS



Basic Information

Catalog Number:

60999-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG23433

GenBank Accession Number:

BC035974

GeneID (NCBI):

1665

UNIPROT ID:

O43143

Full Name:

DEAH (Asp-Glu-Ala-His) box
polypeptide 15

Calculated MW:

795 aa, 91 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F7C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

