

For Research Use Only

FBXL4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61001-1-PBS



Basic Information

Catalog Number: 61001-1-PBS	GenBank Accession Number: BC055010	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 26235	CloneNo.: 2F3F5
Isotype: IgG1	UNIPROT ID: Q9UKA2	
Immunogen Catalog Number: AG28326	Full Name: F-box and leucine-rich repeat protein 4	
	Calculated MW: 621 aa, 70 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51441-1, FBXL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61001-1-PBS. Detection antibody: 61001-2-PBS.



Cytometric bead array standard curve of MP51441-1, FBXL4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61001-1-PBS. Detection antibody: 61001-2-PBS. Standard: Ag28326. Range: 0.781-100 ng/mL