

For Research Use Only

DBP Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61003-1-PBS



Basic Information

Catalog Number:

61003-1-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG14692

GenBank Accession Number:

BC011965

GeneID (NCBI):

1628

UNIPROT ID:

Q10586

Full Name:

D site of albumin promoter (albumin
D-box) binding protein

Calculated MW:

325 aa, 34 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1B2D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51447-1, DBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61003-1-PBS. Detection antibody: 61003-2-PBS.

Cytometric bead array standard curve of MP51447-1, DBP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61003-1-PBS. Detection antibody: 61003-2-PBS. Standard: Ag14692. Range: 0.195-50 ng/mL.