

For Research Use Only

CMIP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61004-2-PBS



Basic Information

Catalog Number: 61004-2-PBS	GenBank Accession Number: BC038113	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 80790	CloneNo.: 2B12A9
Isotype: IgG1	UNIPROT ID: Q8IY22	
Immunogen Catalog Number: AG17875	Full Name: c-Maf-inducing protein	
	Calculated MW: 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51453-1, CMIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61004-1-PBS. Detection antibody: 61004-2-PBS.



Cytometric bead array standard curve of MP51453-1, CMIP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61004-1-PBS. Detection antibody: 61004-2-PBS. Standard: Ag17875. Range: 1.563-400 ng/mL



Bi-layer interferometry (BLI) kinetic assays of 61004-2-Ig against Human CMIP were performed. The affinity constant is below 1 pM.