

For Research Use Only

# ZSCAN4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61006-1-PBS



## Basic Information

Catalog Number:

61006-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30139

GenBank Accession Number:

BC101738

GeneID (NCBI):

201516

UNIPROT ID:

Q8NAM6

Full Name:

zinc finger and SCAN domain  
containing 4

Calculated MW:

433 aa, 49 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D4A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51455-1, ZSCAN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61006-1-PBS. Detection antibody: 61006-2-PBS.



Cytometric bead array standard curve of MP51455-1, ZSCAN4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61006-1-PBS. Detection antibody: 61006-2-PBS. Standard: Ag30139. Range: 0.195-200 ng/mL