

For Research Use Only

CBLL1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61007-2-PBS



Basic Information

Catalog Number:

61007-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30913

GenBank Accession Number:

BC027460

GeneID (NCBI):

79872

UNIPROT ID:

Q75N03

Full Name:

Cas-Br-M (murine) ecotropic retroviral transforming sequence-like 1

Calculated MW:

491 aa, 55 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A1F7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51461-1, CBL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61007-1-PBS. Detection antibody: 61007-2-PBS.



Cytometric bead array standard curve of MP51461-1, CBL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61007-1-PBS. Detection antibody: 61007-2-PBS. Standard:Ag30913. Range: 3.125-400 ng/mL.