

For Research Use Only

BTN3A1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61008-2-PBS



Basic Information

Catalog Number: 61008-2-PBS	GenBank Accession Number: BC121800	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 11119	CloneNo.: 1C12A9
Isotype: IgG1	UNIPROT ID: O00481	
Immunogen Catalog Number: EG1915	Full Name: butyrophilin, subfamily 3, member A1	
	Calculated MW: 513 aa, 58 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51463-1, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-2-PBS.



Cytometric bead array standard curve of MP51463-1, BTN3A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61008-1-PBS. Detection antibody: 61008-2-PBS. Standard: Eg1915. Range: 0.391-50 ng/mL.