

For Research Use Only

Cyclin E2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61016-1-PBS



Basic Information

Catalog Number:

61016-1-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG29625

GenBank Accession Number:

BC020729

GeneID (NCBI):

9134

UNIPROT ID:

O96020

Full Name:

cyclin E2

Calculated MW:

374 aa, 44 kDa

Purification Method:

Protein A purification

CloneNo.:

1E7B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51494-1, CCNE2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61015-1-PBS. Detection antibody: 61016-1-PBS.



Cytometric bead array standard curve of MP51494-1, CCNE2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61015-1-PBS. Detection antibody: 61016-1-PBS. Standard:Ag23601. Range: 6.250-200 ng/mL