

For Research Use Only

TNFSF11/RANKL Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61019-1-PBS



Basic Information

Catalog Number: 61019-1-PBS	GenBank Accession Number: NM_003701	Purification Method: Protein G Magarose purification
Source: Mouse	GeneID (NCBI): 8600	CloneNo.: 1G2A2
Isotype: IgG1	UNIPROT ID: O14788	
Immunogen Catalog Number: EG0982	Full Name: tumor necrosis factor (ligand) superfamily, member 11	
	Calculated MW: 28KD	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP50526-2, RANKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61019-1-PBS. Detection antibody: 61019-2-PBS.



Cytometric bead array standard curve of MP50526-2, RANKL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61019-1-PBS. Detection antibody: 61019-2-PBS. Standard: Eg0982. Range: 1.563-200 ng/mL