

For Research Use Only

RNF2 Monoclonal antibody, PBS Only

Catalog Number: 61023-2-PBS



Basic Information

Catalog Number:

61023-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9129

GenBank Accession Number:

BC012583

GeneID (NCBI):

6045

UNIPROT ID:

Q99496

Full Name:

ring finger protein 2

Calculated MW:

336 aa, 38 kDa

Observed MW:

38 kDa

Purification Method:

Protein G purification

CloneNo.:

3D3B9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

RNF2 (also known as ding, Ring1B or Ring2) is a member of the Ring finger protein family, which functions as E3 ubiquitin ligase for monoubiquitination of histone H2A at lysine 119 (H2AK119ub). Many studies have demonstrated that overexpressed RNF2 was involved in the pathological progression of multiple cancers and has an impact on their clinical features. For instance, the upregulated expression level of RNF2 is positively correlated with the occurrence and progression of hepatocellular carcinoma, melanoma, prostate cancer, breast cancer, pancreatic cancer, gastric cancer, and bladder urothelial carcinoma, as well as with the radioresistance of lung cancer and chemoresistance of ovarian cancer (PMID: 33400276). The molecular weight of RNF2 is 38 kDa.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 61023-2-Ig (RNF2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 61023-2-PBS in a different storage buffer formulation.