

For Research Use Only

# SERPINB4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61039-2-PBS



## Basic Information

<b>Catalog Number:</b> 61039-2-PBS	<b>GenBank Accession Number:</b> BC017401	<b>Purification Method:</b> Protein G Magarose purification
<b>Source:</b> Mouse	<b>GeneID (NCBI):</b> 6318	<b>CloneNo.:</b> 5A2C11
<b>Isotype:</b> IgG1	<b>UNIPROT ID:</b> P48594	
<b>Immunogen Catalog Number:</b> AG1964	<b>Full Name:</b> serpin peptidase inhibitor, clade B (ovalbumin), member 4	
	<b>Calculated MW:</b> 390 aa, 45 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

