

For Research Use Only

CCR9 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61042-1-PBS



Basic Information

Catalog Number: 61042-1-PBS	GenBank Accession Number: BC069678	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 10803	CloneNo.: 1E5H1
Isotype: IgG1	UNIPROT ID: P51686	
Immunogen Catalog Number: AG30128	Full Name: chemokine (C-C motif) receptor 9	
	Calculated MW: 369 aa, 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP5155-1, CCR9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61042-1-PBS. Detection antibody: 61042-2-PBS.



Cytometric bead array standard curve of MP5155-1, CCR9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61042-1-PBS. Detection antibody: 61042-2-PBS. Standard: Ag30128. Range: 1.563-200 ng/mL.