

For Research Use Only

CCR9 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61042-2-PBS



Basic Information

Catalog Number:

61042-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30128

GenBank Accession Number:

BC069678

GeneID (NCBI):

10803

UNIPROT ID:

P51686

Full Name:

chemokine (C-C motif) receptor 9

Calculated MW:

369 aa, 42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B4G1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP5155-1, CCR9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61042-1-PBS. Detection antibody: 61042-2-PBS.



Cytometric bead array standard curve of MP5155-1, CCR9 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61042-1-PBS. Detection antibody: 61042-2-PBS. Standard: Ag30128. Range: 1.563-200 ng/mL.