

For Research Use Only

# Uteroglobulin Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61047-1-PBS



## Basic Information

Catalog Number:

61047-1-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG26315

GenBank Accession Number:

BC004481

GeneID (NCBI):

7356

UNIPROT ID:

P11684

Full Name:

secretoglobin, family 1A, member 1  
(uteroglobulin)

Calculated MW:

10 kDa

Purification Method:

Protein A purification

CloneNo.:

2C10B2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51593-1, Uteroglobin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61047-1-PBS. Detection antibody: 61047-2-PBS.



Cytometric bead array standard curve of MP51593-1, Uteroglobin Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61047-1-PBS. Detection antibody: 61047-2-PBS. Standard: Ag26315. Range: 0.156-20 ng/mL