

For Research Use Only

SMURF2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 61048-1-PBS

Basic Information

Catalog Number:

61048-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25734

GenBank Accession Number:

BC093876

GeneID (NCBI):

64750

UNIPROT ID:

Q9HAU4

Full Name:

SMAD specific E3 ubiquitin protein
ligase 2

Calculated MW:

748 aa, 86 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B12B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS.



Cytometric bead array standard curve of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS. Standard: Ag25734. Range: 0.781-100 ng/mL