

For Research Use Only

SMURF2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 61048-2-PBS

Basic Information

Catalog Number: 61048-2-PBS	GenBank Accession Number: BC093876	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 64750	CloneNo.: 1F6D1
Isotype: IgG1	UNIPROT ID: Q9HAU4	
Immunogen Catalog Number: AG25734	Full Name: SMAD specific E3 ubiquitin protein ligase 2	
	Calculated MW: 748 aa, 86 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS.



Cytometric bead array standard curve of MP51609-1, SMURF2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61048-1-PBS. Detection antibody: 61048-2-PBS. Standard: Ag25734. Range: 0.781-100 ng/mL