

For Research Use Only

RPA3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61050-1-PBS



Basic Information

Catalog Number:

61050-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27853

GenBank Accession Number:

BC005264

GeneID (NCBI):

6119

UNIPROT ID:

P35244

Full Name:

replication protein A3, 14kDa

Calculated MW:

13 kDa

Purification Method:

Protein G purification

CloneNo.:

1D6H11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51618-1, RPA3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61050-1-PBS. Detection antibody: 61050-2-PBS.



Cytometric bead array standard curve of MP51618-1, RPA3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61050-1-PBS. Detection antibody: 61050-2-PBS. Standard: Ag27853. Range: 3.125-400 ng/mL