

For Research Use Only

TMSB4X Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61053-1-PBS



Basic Information

Catalog Number: 61053-1-PBS	GenBank Accession Number: BC001631	Purification Method: Protein A Magarose purification
Source: Mouse	GeneID (NCBI): 7114	CloneNo.: 3D3G1
Isotype: IgG3	UNIPROT ID: P62328	
Immunogen Catalog Number: AG14097	Full Name: thymosin beta 4, X-linked	
	Calculated MW: 63 aa, 7 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

