

For Research Use Only

GGH Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61055-1-PBS



Basic Information

Catalog Number:

61055-1-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28878

GenBank Accession Number:

BC025025

GeneID (NCBI):

8836

UNIPROT ID:

Q92820

Full Name:

gamma-glutamyl hydrolase
(conjugase, foyl/polygammaglutamyl
hydrolase)

Calculated MW:

318 aa, 36 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1E8E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51645-1, GGH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61055-2-PBS. Detection antibody: 61055-1-PBS.



Cytometric bead array standard curve of MP51645-1, GGH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61055-2-PBS. Detection antibody: 61055-1-PBS. Standard: Ag28878. Range: 1.563-200 ng/mL.