

For Research Use Only

GGH Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61055-2-PBS



Basic Information

Catalog Number: 61055-2-PBS	GenBank Accession Number: BC025025	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 8836	CloneNo.: 1A9E3
Isotype: IgG1	UNIPROT ID: Q92820	
Immunogen Catalog Number: AG28878	Full Name: gamma-glutamyl hydrolase (conjugase, foyl/polygammaglutamyl hydrolase)	
	Calculated MW: 318 aa, 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51645-1, GGH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61055-2-PBS. Detection antibody: 61055-1-PBS.



Cytometric bead array standard curve of MP51645-1, GGH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61055-2-PBS. Detection antibody: 61055-1-PBS. Standard: Ag28878. Range: 1.563-200 ng/mL.