

For Research Use Only

NAGLU Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61059-1-PBS



Basic Information

Catalog Number: 61059-1-PBS	GenBank Accession Number: BC053991	Purification Method: Protein A purification
Source: Mouse	GeneID (NCBI): 4669	CloneNo.: 1F1H4
Isotype: IgG3	UNIPROT ID: P54802	
Immunogen Catalog Number: AG10061	Full Name: N-acetylglucosaminidase, alpha-	
	Calculated MW: 743 aa, 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51629-1, NAGLU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61059-1-PBS. Detection antibody: 61059-2-PBS.



Cytometric bead array standard curve of MP51629-1, NAGLU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61059-1-PBS. Detection antibody: 61059-2-PBS. Standard: Ag10061. Range: 1.563-200 ng/mL.