

For Research Use Only

HIF3A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61060-1-PBS



Basic Information

Catalog Number:

61060-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26574

GenBank Accession Number:

NM_022462

GeneID (NCBI):

64344

UNIPROT ID:

Q9Y2N7

Full Name:

hypoxia inducible factor 3, alpha subunit

Calculated MW:

72 kDa

Purification Method:

Protein G purification

CloneNo.:

3A5C6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51633-1, HIF3A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61060-1-PBS. Detection antibody: 61061-1-PBS.



Cytometric bead array standard curve of MP51633-1, HIF3A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61060-1-PBS. Detection antibody: 61061-1-PBS. Standard: Ag26574. Range: 0.781-100 ng/mL.