

For Research Use Only

# AKAP7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61062-2-PBS



## Basic Information

Catalog Number:

61062-2-PBS

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG28790

GenBank Accession Number:

BC016927

GeneID (NCBI):

9465

UNIPROT ID:

O43687

Full Name:

A kinase (PRKA) anchor protein 7

Calculated MW:

81 aa, 9 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E2F2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array sample test of MP51634-1, AKAP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61062-1-PBS. Detection antibody: 61062-2-PBS.



Cytometric bead array standard curve of MP51634-1, AKAP7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61062-1-PBS. Detection antibody: 61062-2-PBS. Standard: Ag28790. Range: 1.563-200 ng/mL.