

For Research Use Only

PFKFB1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61064-1-PBS



Basic Information

Catalog Number: 61064-1-PBS	GenBank Accession Number: BC096079	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 5207	CloneNo.: 1A5E2
Isotype: IgG1	UNIPROT ID: P16118	
Immunogen Catalog Number: AG16752	Full Name: 6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 1	
	Calculated MW: 471 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51639-1, PFKFB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61064-1-PBS. Detection antibody: 61064-2-PBS.



Cytometric bead array standard curve of MP51639-1, PFKFB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61064-1-PBS. Detection antibody: 61064-2-PBS. Standard: Ag16752. Range: 3.125-400 ng/mL.