

For Research Use Only

PTPN6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61065-1-PBS



Basic Information

Catalog Number:

61065-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG21415

GenBank Accession Number:

BC002523

GeneID (NCBI):

5777

UNIPROT ID:

P29350

Full Name:

protein tyrosine phosphatase, non-receptor type 6

Calculated MW:

597 aa, 68 kDa

Purification Method:

Protein G purification

CloneNo.:

1D8E5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51641-1, PTPN6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61065-1-PBS. Detection antibody: 61065-2-PBS.



Cytometric bead array standard curve of MP51641-1, PTPN6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61065-1-PBS. Detection antibody: 61065-2-PBS. Standard: Ag21415. Range: 1.563-200 ng/mL.