

For Research Use Only

MAN2B2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61069-2-PBS



Basic Information

Catalog Number:

61069-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29258

GenBank Accession Number:

BC033307

GeneID (NCBI):

23324

UNIPROT ID:

Q9Y2E5

Full Name:

mannosidase, alpha, class 2B, member 2

Calculated MW:

1009 aa, 114 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A10C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51660-1, MAN2B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61069-1-PBS. Detection antibody: 61069-2-PBS.



Cytometric bead array standard curve of MP51660-1, MAN2B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61069-1-PBS. Detection antibody: 61069-2-PBS. Standard: Ag29258. Range: 0.781-100 ng/mL