

For Research Use Only

# WNT2B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61070-1-PBS



## Basic Information

Catalog Number:

61070-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28705

GenBank Accession Number:

BC141825

GeneID (NCBI):

7482

UNIPROT ID:

Q93097

Full Name:

wingless-type MMTV integration site  
family, member 2B

Calculated MW:

391 aa, 44 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F11C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51661-1, WNT2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61070-1-PBS. Detection antibody: 61070-2-PBS.



Cytometric bead array standard curve of MP51661-1, WNT2B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61070-1-PBS. Detection antibody: 61070-2-PBS. Standard: Ag28705. Range: 3.125-400 ng/mL