

For Research Use Only

# SOX6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61071-1-PBS



## Basic Information

Catalog Number:

61071-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5794

GenBank Accession Number:

BC047064

GeneID (NCBI):

55553

UNIPROT ID:

P35712

Full Name:

SRY (sex determining region Y)-box 6

Calculated MW:

92 kDa

Purification Method:

Protein G purification

CloneNo.:

1H5A6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:



T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array sample test of MP51662-1, SOX6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61071-1-PBS. Detection antibody: 61071-2-PBS.

Cytometric bead array standard curve of MP51662-1, SOX6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61071-1-PBS. Detection antibody: 61071-2-PBS. Standard: Ag5794. Range: 0.781-100 ng/mL