

For Research Use Only

SOX6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61071-2-PBS



Basic Information

| | | |
|--|---|--|
| Catalog Number: 61071-2-PBS | GenBank Accession Number: BC047064 | Purification Method: Protein G Magarose purification |
| Source: Mouse | GeneID (NCBI): 55553 | CloneNo.: 2B12C1 |
| Isotype: IgG1 | UNIPROT ID: P35712 | |
| Immunogen Catalog Number: AG5794 | Full Name: SRY (sex determining region Y)-box 6 | |
| | Calculated MW: 92 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51662-1, SOX6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61071-1-PBS. Detection antibody: 61071-2-PBS.



Cytometric bead array standard curve of MP51662-1, SOX6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61071-1-PBS. Detection antibody: 61071-2-PBS. Standard: Ag5794. Range: 0.781-100 ng/mL