

For Research Use Only

HMGN2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61074-1-PBS



Basic Information

Catalog Number:

61074-1-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28531

GenBank Accession Number:

BC014644

GeneID (NCBI):

3151

UNIPROT ID:

P05204

Full Name:

high-mobility group nucleosomal binding domain 2

Calculated MW:

9 aa, 3 kDa

Purification Method:

Protein G purification

CloneNo.:

1B9D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51668-1, HMGN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS.



Cytometric bead array standard curve of MP51668-1, HMGN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS. Standard:Ag28531. Range: 0.781-100 ng/mL



Biolayer interferometry (BLI) kinetic assays of 61074-1-Ig against Human HMGN2 were performed. The affinity constant is 1.87 nM.