

For Research Use Only

HMGN2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61074-2-PBS



Basic Information

Catalog Number: 61074-2-PBS	GenBank Accession Number: BC014644	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 3151	CloneNo.: 2D5D3
Isotype: IgG1	UNIPROT ID: P05204	
Immunogen Catalog Number: AG28531	Full Name: high-mobility group nucleosomal binding domain 2	
	Calculated MW: 9 aa, 3 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array sample test of MP51668-1, HMGN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS.



Cytometric bead array standard curve of MP51668-1, HMGN2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61074-1-PBS. Detection antibody: 61074-2-PBS. Standard:Ag28531. Range: 0.781-100 ng/mL



Bi-layer interferometry (BLI) kinetic assays of 61074-2-Ig against Human HMGN2 were performed. The affinity constant is 1.34 nM.