

For Research Use Only

# TBX3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61079-1-PBS



## Basic Information

Catalog Number: 61079-1-PBS	GenBank Accession Number: BC025258	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 6926	CloneNo.: 2B2A9
Isotype: IgG1	UNIPROT ID: O15119	
Immunogen Catalog Number: AG10195	Full Name: T-box 3	
	Calculated MW: 79 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51688-1, TBX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61079-1-PBS. Detection antibody: 61079-2-PBS.



Cytometric bead array standard curve of MP51688-1, TBX3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61079-1-PBS. Detection antibody: 61079-2-PBS. Standard: Ag10195. Range: 0.781-100 ng/mL.