

For Research Use Only

CTAK1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61082-3-PBS



Basic Information

Catalog Number:

61082-3-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8303

GenBank Accession Number:

BC024773

GeneID (NCBI):

4140

UNIPROT ID:

P27448

Full Name:

MAP/microtubule affinity-regulating kinase 3

Calculated MW:

78 kDa, 80 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F12B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51716-2, CTAK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61082-1-PBS. Detection antibody: 61082-3-PBS.



Cytometric bead array standard curve of MP51716-2, CTAK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61082-1-PBS. Detection antibody: 61082-3-PBS. Standard: Ag8303. Range: 3.125-400 ng/mL