

For Research Use Only

# SNAPC2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61084-2-PBS



## Basic Information

Catalog Number:

61084-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9689

GenBank Accession Number:

BC011868

GeneID (NCBI):

6618

UNIPROT ID:

Q13487

Full Name:

small nuclear RNA activating complex, polypeptide 2, 45kDa

Calculated MW:

334 aa, 36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G8B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array sample test of MP51708-1, SNAPC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61084-1-PBS. Detection antibody: 61084-2-PBS.



Cytometric bead array standard curve of MP51708-1, SNAPC2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61084-1-PBS. Detection antibody: 61084-2-PBS. Standard: Ag9689. Range: 3.125-400 ng/mL