

For Research Use Only

ADORA2A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61087-2-PBS



Basic Information

Catalog Number:

61087-2-PBS

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG11626

GenBank Accession Number:

BC013780

GeneID (NCBI):

135

UNIPROT ID:

P29274

Full Name:

adenosine A2a receptor

Calculated MW:

412 aa, 45 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2G8C1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51717-1, ADORA2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61087-1-PBS. Detection antibody: 61087-2-PBS.



Cytometric bead array standard curve of MP51717-1, ADORA2A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61087-1-PBS. Detection antibody: 61087-2-PBS. Standard:Ag11626. Range: 3.125-400 ng/mL.