

For Research Use Only

EML4 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 61094-2-PBS



Basic Information

Catalog Number: 61094-2-PBS	GenBank Accession Number: BC140845	Purification Method: Thiophilic affinity chromatograph
Source: Mouse	GeneID (NCBI): 27436	CloneNo.: 2D6H1
Isotype: IgM	UNIPROT ID: Q9HC35	
Immunogen Catalog Number: AG18755	Full Name: echinoderm microtubule associated protein like 4	
	Calculated MW: 981 aa, 109 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51727-1, EML4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61094-1-PBS. Detection antibody: 61094-2-PBS.



Cytometric bead array standard curve of MP51727-1, EML4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61094-1-PBS. Detection antibody: 61094-2-PBS. Standard: Ag18755. Range: 0.781-100 ng/mL.