

For Research Use Only

PRKDC Monoclonal antibody, PBS Only (Capture)

Catalog Number: 61098-1-PBS



Basic Information

Catalog Number: 61098-1-PBS	GenBank Accession Number: NM_006904	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 5591	CloneNo.: 2E4D5
Isotype: IgG1	UNIPROT ID: P78527	
Immunogen Catalog Number: AG29173	Full Name: protein kinase, DNA-activated, catalytic polypeptide	
	Calculated MW: 469 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Cytometric bead array sample test of MP51771-1, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-3-PBS.



Cytometric bead array sample test of MP51771-2, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-2-PBS.



Cytometric bead array standard curve of MP51771-1, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-3-PBS. Standard:Ag29173. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP51771-2, PRKDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 61098-1-PBS. Detection antibody: 61098-2-PBS. Standard:Ag29173. Range: 0.156-20 ng/mL